Thickness: 4 mil (0,10 mm)
Length: 11.5" (295 mm)

SAME GLOVE - BETTER FOR THE PLANET ChemRest® Series EBT by SHOWA®:
The environmental impact is the only difference SHOWA® gloves with Eco Best Technology® (EBT) are engineered for accelerated biodegradation in biologically active landfills. Independent certified laboratories performing long-term landfill biodegradation testing according to ASTM D5526-12 reported that SHOWA's gloves with EBT achieved 82.0% biodegradation in only 386 days, while gloves without EBT achieved only 1.9% biodegradation over the same period of time. These results may not be indicative of future biodegradation.

POWDER-FREE • BIODEGRADABLE • SINGLE USE • NITRILE GLOVES

EN ISO 21420:2020
Level 5 (Min. 1 - Max. 5)

Dexterity	Handvaardigheid	Kätevyyys
Dextérité	Böjlighet	Sprawność manualna
Fingerfertigkeit	Bevægelighed	Yetenekler
Destreza	Bevegelighet	البراعة اليدوية

EN ISO 374-1:2016
+A1:2018/Type B

Chemical risks	Risques chimiques	Kemiska risker
Risques chimiques	Chemische Risiken	Kemiske risici
Rischi chimici	Riesgos químicos	Kjemiske risikoer
Riesgos químicos	Riscos Químicos	Kemialliset riskit
Riscos Químicos	Chemische risico's	Zagrożenia chemiczne
		Kimyasal Riskler
		المخاطر الكيميائية

EN 16523-1: 2015+A1:2018 Permeation	EN ISO 374-4: 2019 Degradation
Level 4	24.3%
Level 6	-8.1%
Level 2	29.3%
Level 4	25.6%
Level 6	21.7%

(Min. 1 - Max. 6)

J N-Heptane (142-82-5)	Level 4
K Sodium hydroxide 40% (1310-73-2)	Level 6
O Ammonium hydroxide 25% (1336-21-6)	Level 2
P Hydrogen peroxide 30% (7722-84-1)	Level 4
T Formaldehyde 37% (50-00-0)	Level 6

Level 1 > 10 min	Level 3 > 60 min	Level 5 > 240 min
Level 2 > 30 min	Level 4 > 120 min	Level 6 > 480 min

EN ISO 374-5:2016

Micro-organisms risks	Risques liés aux micro-organismes	Risiker för mikroorganismer
Risques liés aux micro-organismes	Risiko durch Mikroorganismen	Risici for mikroorganismer
Risiko durch Mikroorganismen	Rischi di microorganismi	Risikoer for mikroorganismer
Rischi di microorganismi	Riesgos de microorganismos	Mikro-organismien riskit
Riesgos de microorganismos	Riscos de micro-organismos	Ryzyko drobnoustrojów
Riscos de micro-organismos	Risiko's van micro-organismen	Mikroorganizma riskleri
Risiko's van micro-organismen		مخاطر الكائنات المجهرية

Protection against bacteria and fungi - Pass
Protection against viruses - Pass

EN ISO 374-2:2019
Penetration

Level	AQL	Inspection level
Level 3	<0.65	G1
Level 2	<1.5	G1
Level 1	<4.0	S4

For Contact with Foodstuffs: COMMISSION REGULATION (EU) No. 10/2011, COMMISSION REGULATIONS (EC) No. 1935/2004, No. 2023/2006 and FDA 21 CFR 177.2600

CE 2777

Regulation (EU) 2016/425

Module B EU Type Examination and Module C2
Ongoing Conformity SATRA Technology Europe Ltd.
Bracetown Business Park, Clonee, Dublin D15 YN2P,
Ireland (Notified Body No. 2777)

EN Gloves have a shelf life of 5 years • Potential allergen: Carbamate • Gloves provide protection from chemical hazards shown • Before usage, inspect the gloves for any defects or imperfections • It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion, and degradation • The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only (except in cases where the glove is equal to or over 400 mm - where the cuff is tested also) and relates only to the chemical tested, it can be different if the chemical is used in a mixture • Permeation information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals • When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties • Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly • Where the test specimens gave an increased puncture force after chemical exposure, the result is reported as a negative degradation • For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves • The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen • Donning: Thoroughly wash hands. Select the appropriately sized gloves. Hold with one hand and insert the other. Pull glove cuff towards wrist to cover as much skin as possible and secure glove. Check to make sure there are no holes or tears. Doffing: Grasp the outside edge of the glove near the wrist. Peel the remaining glove off from the inside, creating a "bag" containing both gloves. Peel the glove away from the hand, turning it inside out. • All Showa products must be stored unopened at room temperature and protected from heat above 40°C (104°F), humidity, sunlight, ozone, pests and sharp objects • It is recommended that all Showa products be subject to a "First In, First Out" stock rotation. Products with or without an expiry date do not release the customer from ensuring product viability upon receipt and do not warrant the fitness of a product for any particular use • Discard used gloves in compliance with local regulations

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<p>EU DECLARATION OF CONFORMITY</p> <p>www.showagroup.com</p>	 <p>ChemRest[®]</p> <p>Chemical Resistant Glove Directory</p> <p>KNOW YOU'RE PROTECTED</p>	<p>www.chemrest.com</p> <p>Type on the icon is in accordance with EN ISO 374-1:2016+A1:2018 standard.</p>
<p> Manufacturer</p> <p>Showa Best Glove, Inc.</p> <p>579 Edison Street</p> <p>Menlo, GA 30071 USA</p> <p>706.862.2302</p> <p>800.241.0323</p> <p>USA@SHOWAgroup.com</p>	<p> Importer & Authorised Representative</p> <p>Showa International (Netherlands) B.V.</p> <p>WTC Amsterdam, Tower 2, 3rd floor</p> <p>Strawinskylaan 1817, 1077 XX Amsterdam</p> <p>The Netherlands</p> <p>+31 (0)88-004 2100</p> <p>info@SHOWAgroup.eu</p>	