SHOWA
ChemRest
5500PF

6/XS • 6-7/S • 7-8/M • 8-9/L • 9-10/XL • 10-11/XXL

POWDER-FREE • NON-STERILE • SINGLE USE • NITRILE GLOVES

Thickness: 4 mil(0.10 mm)

length: 9.5" (240 mm)

Regulation (EU) 2016/425

Module B EU Type Examination and Module D Ongoing Conformity

SATRA Technology Europe Ltd.

Bracetown Business Park, Clonee, Meath D15 YN2P,

Ireland (Notified Body No. 2777)

For Contact with Foodstuffs
COMMISSION REGULATION (EC) No.1935/2004

EU DECLARATION OF CONFORMITY
www.showagroup.com

-Authorized Representative / Importer

EU SHOWA International (Netherlands) B.V.

WTC Tower I, Strawinskylaan 1817,
1077 XX Amsterdam, The Netherlands

+31 (0)88-004 2100

info@SHOWAgroup.eu

-Manufacturer

JP SHOWA GLOVE Co.

565 Tohori, Himeji, Hyogo, 670-0802 Japan

Made in China

TO FULLY ASSESS YOUR ASSOCIATED RISKS
WITH CHEMICAL SUBSTANCES USE :

www.chemrest.com

Type on the icon in accordance with
EN ISO 374-1:2016+A1:2018 standard.

[EN]

• These gloves have an obsolescence of 3 years with correct storage in an unused condition. It cannot be guaranteed as the material will be affected by external factors such as temperature, humidity, ozone, etc. • Potential allergen: Carbamate. • Gloves provide protection from chemical hazards shown. • Before usage, inspect the gloves for any defects or imperfections. • It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion, and degradation. • The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only (except in cases where the glove is equal to or over 400 mm - where the cuff is tested also) and relates only to the chemical tested, it can be different if the chemical is used in a mixture. • Permeation information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals. • When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. • Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. • Where the test specimens gave an increased puncture force after chemical exposure, the result is reported as a negative degradation. • For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves. • The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen. • Donning: Thoroughly wash hands. Select the appropriately sized gloves. Hold with one hand and Insert the other. Pull glove cuff towards wrist to cover as much skin as possible and secure glove. Check to make sure there are no holes or tears. Doffing: Grasp the outside edge of the glove near the wrist. Peel the remaining glove off from the inside, creating a "bag" containing both gloves. Peel the glove away from the hand, turning it inside out. • All Showa products must be stored unopened at room temperature and protected from heat above 40°C (104°F), humidity, sunlight, ozone, pests and sharp objects. • It is recommended that all Showa products be subject to a "First In, First Out" stock rotation. Products with or without an expiry date do not release the customer from ensuring product viability upon receipt and do not warrant the fitness of a product for any particular use. • Discard used gloves in compliance with local regulations.

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38.3%

24.7%

Degradation levels indicate the change in puncture resistance of the gloves after exposure to the challenge chemical.

EN ISO 374-5:2016

EN ISO 374-5

VIRUS

Micro-organisms risks
Risques liés aux micro-organismes
Risiko durch Mikroorganismen
Rischi di microrganismi
Riesgos de microorganismos
Riscos de micro-organismos
Risiko's van micro-organismen

Risker för mikroorganismer
Risici for mikroorganismer
Risikoer for mikroorganismer
Mikro-organismien riskit
Ryzyko drobnoustrojów
Mikroorganizma riskleri
مخاطر الكائنات المجهرية

Protection against bacteria and fungi - Pass

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