SHOWA
ChemRest
7665PF ESD

6/XS • 7/S • 8/M • 9/L • 10/XL • 11/XXL

POWDER-FREE • NON-STERILE • SILICONE-FREE • SINGLE USE • ANTI STATIC • NITRILE GLOVES

Silicone Free refers; Silicone containing materials are not intentionally used in the manufacturing process of this product.

Thickness: 6 mil(0.15 mm)

Length: 12" (300 mm)

Regulation (EU) 2016/425

CE 0598
CAT.III

Module B EU Type Examination Certificate, Notified Body 2777
SATRA Technology Europe Ltd.
Bracetown Business Park, Clonee, Meath, D15 YN2P, Ireland

Module D Ongoing Conformity, Notified Body 0598
SGS Fimko Oy
Takomotie 8, FI-00380 Helsinki, Finland

EU DECLARATION OF CONFORMITY
www.showagroup.com

-Authorized Representative / Importer

EU SHOWA International (Netherlands) B.V.
WTC Tower I, Strawinskylaan 1817,
1077 XX Amsterdam, The Netherlands
+31 (0)88-004 2100
info@SHOWAgroup.eu

-Manufacturer

US SHOWA Best Glove Inc.
579 Edison Street, Menlo, GA 30731 USA

Made in Sri Lanka

TO FULLY ASSESS YOUR ASSOCIATED RISKS
WITH CHEMICAL SUBSTANCES USE :
www.chemrest.com
Type on the icon is in accordance with
EN ISO 374-1 standard.

[EN]
• These gloves have an obsolescence of 5 years with correct storage in an unused condition. It cannot be guaranteed as the material will be affected by external factors such as temperature, humidity, ozone, etc. • In the case of allergic reaction, medical aid should be sought immediately. • For single use only. • Protection against chemical risks, with electrostatic properties. • Do not use where there are electrical, thermal or entanglement risks. • Gloves provide protection from chemical hazards shown. •Before usage, inspect the gloves for any defects or imperfections. • It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation. • When used, protective gloves may provide less resistance to the dangerous chemicals due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistance gloves. • Permeation information does not reflect the actual duration of protection in the workplace due to other factors that may influence the performance and the differentiation between mixtures and pure chemicals. • The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only and relates only to the chemical tested. It can be different if the chemical is used in a mixture. • Designed to protect against micro-organisms and comply with EN ISO 374-5 requirements. • The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen. • Not tested against viruses. • When putting on gloves, grasp the glove at the bottom of the cuff and pull the glove over your hand. Ensure the gloves fit well. Repeat on the other hand. • When removing gloves, place the fingertips into the palm of the other glove. Pull the glove until off. Repeat on the other hand. • All Showa products must be stored unopened at room temperature and protected from heat above 40° C (104° F), humidity, sunlight, ozone, pests and sharp objects. • Discard used gloves in compliance with local regulations. [Warning] • The person wearing the electrostatic dissipative protective gloves shall be properly earthed e.g. by wearing adequate footwear. • Electrostatic dissipative protective gloves shall not be unpacked, opened, adjusted or removed whilst in flammable or explosive atmospheres or while handling flammable or explosive substances. • The electrostatic properties of the protective gloves might be adversely affected by ageing, wear, contamination and damage, and might not be sufficient for oxygen enriched flammable atmospheres where additional assessments are necessary.

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