

SHOWA M7002PF

SHOWA M7002PF is an FDA Class 1 Medical Device intended for medical purposes that is worn on the examiner's hands to prevent contamination between patient and examiner. By using SHOWA's Eco Best Technology®, the M7002PF's are completely eco-friendly, breaking down by 82% in just 386 days when tested in a laboratory.



BENEFITS

- Disposable
- Latex-free
- Increased dexterity
- Flexible
- Textured Fingertips

FEATURES

- Eco Best Technology® (EBT)
- Powder-free
- Rolled Cuff

NORMS & CERTIFICATES

INDUSTRIES



Healthcare



Laboratory



Municipal
Services



Pharmaceutical

HAZARDS



Medical



Single Use

TRADES & APPLICATIONS

- Hospital & medical care
- Emergency services
- Police & defense
- Laboratory, pharma & analysis

GET IN TOUCH

SHOWA International | 579 Edison Street, Menlo, GA 30731, USA | showagroup.com |

USA@SHOWAGroup.com

© SHOWA GROUP 2025 | All rights reserved

PACKAGING

- Dispenser par case: 10
- Gloves per dispenser: 200
- Gloves per case: 2000
- Packaging info - XXL size:
& XL size: 180 gloves per dispenser

THICKNESS

3 mil

LENGTH

9.5"

FDA Class I Medical Device

COATING

- Nitrile
- Biodegradable

LINERS

- Disposable
- Unsupported

SIZES

6/XS | 7/S | 8/M | 9/L | 10/XL |
11/XXL

COLOUR

- Blue

MATERIAL

- Single use
- Unsupported

GRIP

- Textured

USER INSTRUCTIONS

Gloves provide protection from chemical and mechanical hazards shown. Do not use gloves that show signs of wear. If required, cleanse outer surface of glove with running water. Discard used gloves in compliance with local regulations. Do not wear gloves when there is a risk of entanglement by moving parts of machines.

DISCLAIMER

The descriptions, characteristics, applications and photos are given for information purposes and do not constitute a contractual commitment. The manufacturer reserves the right to make any modifications it deems necessary.

GET IN TOUCH

SHOWA International | 579 Edison Street, Menlo, GA 30731, USA | showagroup.com |
USA@SHOWAGroup.com
© SHOWA GROUP 2025 | All rights reserved